

High-throughput Sodium Flux Assay for Na_v1.5 using Assay Ready Cells

Shu Z. Wiley

ION Biosciences, San Marcos, TX 78666, United States

Abstract

Voltage-gated sodium channels (Na_vs) present tremendous therapeutic potential for treating pain, epilepsy, neuromuscular disorders, cardiovascular diseases, respiratory disorders, and cancer. Na_v1.5 is the primary Na_v channel in the heart, and the well-established molecular target for antiarrhythmic drugs and local anesthetics. A fluorescence-based high-throughput sodium flux assay was developed to screen small molecule modulators targeting Na_v1.5. This assay utilizes ION Biosciences' Brilliant Sodium Assay and acCELLerate's assay ready CHO-NaV1.5-DUO™ instaCELLs to provide an effective, simple, and fast method to identify Na_v1.5 activators/inhibitors.

Introduction

Voltage-gated sodium channels (Na_v1.X) are expressed in the cell membranes of excitable and non-excitable cells. In humans there are nine different Na_v α subunits, Na_v1.1 – Na_v1.9. Membrane depolarization facilitates a conformational change of the α subunit, resulting in the opening of the Na⁺-selective channel and inward current of Na⁺. Na_v channels work in concert to initiate action potentials; therefore, it is no surprise that their dysfunction plays a role in numerous neuronal and cardiac disorders. Sodium channel blockers have been developed as anticonvulsant, antiarrhythmic, and local anesthetic drugs in the past. Currently, enormous efforts are underway to find and study the effects of subtype-selective Na_v modulators in pre-clinical disease models [1].

Na_v1.5 channels, encoded by gene SCN5A, are mostly expressed in human heart tissue. Na_v1.5 channel mutations are associated with many cardiac dysfunctions including Brugada syndrome, long QT syndrome, and atrial fibrillation [2]. Na_v1.5 is also a therapeutic target for antiarrhythmic drugs and local anesthetics such as lidocaine and tetracaine [3], and is part of the Comprehensive in vitro Proarrhythmia Assay (CiPA) panel for evaluation of cardiac risk. Thus, there is high interest in screening for Na_v1.5 modulators.

Here we present a fluorescence-based high-throughput sodium flux assay that is suitable for screening for modulators targeting Na_v1.5. This assay utilizes ION Biosciences' Brilliant Sodium Assay and acCELLerate's assay ready CHO-NaV1.5-DUO™ instaCELLs to provide a rapid and simple method to identify Na_v1.5 activators/inhibitors in less than 24 hrs. ION Biosciences Brilliant Sodium Assay contains ION Natrium Green-2 (ING-2), a sodium selective-indicator, and all other reagents necessary for running

a HTS assay. CHO-Nav1.5-DUO™ instaCELLS are qualified frozen aliquots of cells from a recombinant cell line which stably express Nav1.5. The cell line had been developed and validated by B'SYS GmbH, Switzerland and is distributed by acCELLerate in the assay ready instaCELL format. Instantly after thawing and without prior cultivation, CHO-Nav1.5-DUO™ instaCELLS display a high viability (>90%) and exhibit a strong functional expression of Nav1.5.


Experimental overview

Thaw assay ready CHO-Nav1.5-DUO™ instaCELLS and dilute into the desired concentration with complete cell culture media. Seed cells into black wall/clear bottom multi-well plates (30,000 cells/well for 96-well plate and 15,000 cells/well for 384-well plate) and incubate overnight to allow the cells to attach. The next day, make a 1X dye loading solution according to the Brilliant Sodium Assay protocol, remove cell culture media and add dye loading solution into each cell-containing well (100µL/well for 96-well plate and 20µL/well for 384-well plate). Incubate the plates for 1 hr. During the incubation, prepare a 3X concentrated compound plate. Finally, perform assay readout on a kinetic plate reader using the parameters listed below.

Assay format	Instrument	Ex*, Em	Frequency	Baseline read	Drug addition	Post-addition read
96-well	Flexstation 3	515 nm, 545 nm	1 Hz	20 s	50 µL	130 s
384-well	Panoptic	518 nm, 562 nm	1 Hz	20 s	10 µL	300 s

*ING-2 Na⁺ indicator also works with standard 488nm light sources.

When testing for Nav1.5 inhibitors, add inhibitors 15 min prior to your read, then introduce a high concentration of veratridine, a blocker of Nav_v inactivation, during the kinetic read. If a kinetic plate reader isn't available, the assay can be conducted in endpoint mode using a standard plate reader or fluorescence microscope using either GFP or YFP filters.




Materials and reagents

Category	Name	Description	Cat. #	Supplier
Assay reagents	Brilliant Sodium Assay Flex Kit	Assay kit includes Na ⁺ indicator ING-2 and all other reagents needed (e.g., assay buffer, DySolv, Probenecid, and TRS)	9000-10	ION Biosciences
Cell culture	CHO-Nav1.5-DUO™	Assay ready CHO-K1 cells that overexpress Nav1.5	RE305	acCELLerate/B'SYS GmbH
	Advanced MEM	Base media (Complete media contain 5% FBS and 1X GlutaMAX)	12491013	ThermoFisher
	Fetal Bovine Serum (FBS)	Fetal Bovine Serum	35-016-CV	Corning
	GlutaMAX™ Supplement	200 mM L-alanyl-L-glutamine dipeptide in 0.85% NaCl	35050061	ThermoFisher
Tool compounds	Veratridine	Voltage-gated Na ⁺ channel activator		
	Tetracaine	Na ⁺ channel blocker		


Results

Veratridine (VTD) is a steroidal alkaloid neurotoxin isolated from the seeds of the Central America and Mexican plant sabadilla lily. Veratridine causes sodium channels to stay open during a sustained membrane depolarization by abolishing inactivation [4]. Because of its ability to open Na_v channels, veratridine is used in cell-based drug screening assays in which controlling the membrane voltage is




impractical [5]. Fluorescent Na⁺ indicators that measure the intracellular ion concentrations have been successfully used to assess veratridine pharmacological functions on Na_v1.1-Na_v1.7 channels [6]. In our assay, a dose-dependent increase of intracellular fluorescence is detected after veratridine treatment induces an influx of Na⁺ (Fig 2a). Dose-response curve yield an EC₅₀ of 101.8 μM (Fig 2b), which is comparable to previously reported values for veratridine [6,7]. Fluorescence images acquired after the kinetic read also show a sustained elevation in fluorescence intensity at higher veratridine concentrations (Fig 2c). Collectively, these data demonstrate that a sodium flux assay using CHO-Na_v1.5-DUO™ instaCELLs is an accessible and rapid approach to detect Na_v1.5 channel activators.

Tetracaine is an amino-ester class local anesthetic used to numb the eyes, nose, or throat. Tetracaine inhibits Na_v function by blocking the channel at the resting state [8]. To measure tetracaine inhibition, CHO-Na_v1.5-DUO™ instaCELLs were first incubated with different concentrations of tetracaine for 15 min at 37°C. Then 300 μM veratridine was added to wells during kinetic data acquisition. A tetracaine concentration-dependent reduction in fluorescence (F/F₀) was observed in the normalized kinetic data when compared to a veratridine only control (Fig 3a). A tetracaine IC₅₀ of 84.4 μM was calculated (Fig 3b), which is comparable to previously reported values [6]. Fluorescence images acquired after the kinetic read showed cellular fluorescence was lower in the presence of higher concentrations of tetracaine (>333 μM) (Fig 3c). These data demonstrate that CHO-Na_v1.5-DUO™ cell-based sodium flux assay can also be used to identify inhibitors of Na_v1.5 channels.



The 96-well assay was miniaturized to a 384-well format using Wavefront Biosciences Panoptic kinetic plate reader. The assay was tested for screening both Nav_v1.5 channel activators and Nav_v1.5 inhibitors. For activator screening, a quarter of the 384-well plate was treated with 300 μM veratridine as the positive control group and another quarter of the 384-well plate was a buffer only negative control group. Normalized kinetic data shows a large dynamic window of 1.9-fold at 300s after veratridine addition (Fig 3a). The calculated Z' factor of 0.83 indicates a high-quality screening assay to identify Nav_v1.5 channel activators (Fig 3b). For inhibitor screening, a quarter of the 384-well plate was pre-treated with 1 mM tetracaine as the negative control group and another quarter of the 384-well plate was pre-treated with buffer only as the positive control group. High concentration (300 μM) veratridine was added to all wells after the 20s baseline data collection during the kinetic readout. Normalized kinetic data showed a significant dynamic window at 320s timepoint (Fig 3c). The calculated Z' factor of 0.70 indicates a high-quality screening assay to identify Nav_v1.5 channel inhibitors (Fig 3d).



Conclusion

High-throughput sodium flux assays for Na_v1.5 channels using assay ready instaCELLs is a high-quality assay that can be used to identify Na_v1.5 modulators and assess potency of Na_v1.5 activators and inhibitors. Veratridine EC₅₀ and tetracaine IC₅₀ are comparable to reported values, and both activator and inhibitor screens generated Z' factors larger than 0.5. In this assay, veratridine is used instead of membrane depolarization to open sodium channels as part of a high-throughput intracellular Na⁺ concentration screening strategy to identify modulating compounds. The combination of ION Biosciences' Brilliant Sodium Assay and acCELLerate's assay ready cells provides a robust method for testing Na_v1.5 and can be extended to identify other Na_v channel modulators as well.

Acknowledgements

This work was done in collaboration with:



We thank acCELLerate for providing CHO-NaV1.5-DUO™ instaCELLs for this work. CHO-NaV1.5-DUO™ instaCELLs are qualified frozen aliquots of cells from a recombinant cell line which stably express Na_v1.5. The cell line had been developed and validated by B'SYS GmbH, Switzerland and is distributed by acCELLerate in the assay ready instaCELL format.

References

1. Ruiz ML, and Kraus R, Voltage-Gated Sodium Channels: Structure, Function, Pharmacology, and Clinical Indications. *J. Med. Chem.*, **2015**, 24, 7093-7118.
2. Antzelevitch C., Nesterenko V, Shryock JC, Rajamani S, Song Y, and Belardinelli L, The role of late I Na in development of cardiac arrhythmias. *Handb. Exp. Pharmacol.*, **2014**, 221, 137–168.
3. Fozzard HA, Sheets MF, and Hanck DA, The sodium channel as a target for local anesthetic drugs. *Front. Pharmacol.*, **2011**, 2, 68.
4. Ulbricht W, Effects of veratridine on sodium currents and fluxes. *Rev. Physiol. Biochem. Pharmacol.*, **1998**, 133, 1-54.
5. Felix JP, Williams BS, Priest BT, Brochu RM, Dick IE, Warren VA, Yan L, Slaughter RS, Kaczorowski GJ, Smith MM, and Garcia ML, Functional assay of voltage-gated sodium channels using membrane potential sensitive dyes. *Assay Drug Dev. Technol.*, **2004**, 2, 260-268.
6. Tay B, Stewart TA, Davis FM, Deuis JR, and Vetter I, Development of a high-throughput fluorescent no-wash sodium influx assay. *PLoS One.*, **2019**, 14, e0213751.
7. Vetter I, Mozar CA, Durek T, Wingerd JS, Alewood PF, Christie MJ, et al. Characterisation of Na(v) types endogenously expressed in human SH-SY5Y neuroblastoma cells. *Biochem. Pharmacol.*, **2012**, 83, 1562-1571.

8. Li HL, Galue A, Meadows L, and Ragsdale DS, A molecular basis for the different local anesthetic affinities of resting versus open and inactivated states of the sodium channel. *Mol. Pharmacol.*, **1999**, 55, 134-141.
9. Zhang J, Chung TD, and Oldenburg KR, A simple statistical parameter for use in evaluation and validation of high throughput screening assays. *J. Biomol. Screen*, **1999**, 4, 67-73.