

## Brilliant Sodium Assay

Table 1 Kit Contents		<i>Flex Kit</i> Cat# 9000-10		<i>pIONeer Kit</i> Cat# 9000-2		
Label	Name	Size	Qty	Size	Qty	Storage
Reagent A	ING-2 AM	50 µg Vial	10	50 µg Vial	2	-20° C
Reagent B	DMSO <sup>1</sup>	225 µL Vial	1	Not included in pIONeer		-20° C
Reagent C	50X DySolv	4 mL Bottle	1	800 µL Vial	1	4° C
Reagent D	10X Brilliant Assay Buffer	20 mL Bottle	1	4 mL Bottle	1	4° C
Reagent E	50X TRS <sup>2</sup>	4 mL Bottle	1	800 µL Vial	1	4° C
Reagent F	50X Probenecid <sup>2</sup>	4 mL Bottle	1	800 µL Vial	1	4° C

### *Express (100-Plate) Kit - Cat# 9000-100*

Brilliant Sodium Assay is also available in a convenient 100-Plate size kit, perfect for large scale assays and high-throughput screening. These kits are customizable to your institution's needs so that you receive only the reagents you require, in a format that works best for your application. Contact ION Biosciences Sales to discuss how we can tailor our Express kits to fit your needs.

### Description

Sodium (Na<sup>+</sup>) is one of the most important monovalent metal cations in living organisms. Na<sup>+</sup> channels, Na<sup>+</sup>-permeable non-selective monovalent cation channels, and Na<sup>+</sup>-coupled transporters play critical roles including modulating neuronal activity, powering transport of nutrients and signaling molecules, and regulating solute balance. Na<sup>+</sup>-permeable channel and Na<sup>+</sup> transporter-targeted drugs provide effective treatments for a diversity of indications: epilepsy, pain, bipolar disorder, depression, diuresis, and many others. As a result, interest in Na<sup>+</sup>-permeable channels and Na<sup>+</sup> transporters as drug targets remains high.

ION's Brilliant Sodium Assay is a total assay solution for multi-well plate-based, high-throughput measurements of changes in intracellular Na<sup>+</sup> mediated through a wide-variety of plasma membrane and intracellular sodium channels and transporters. In multi-well, plate-based formats, the Brilliant Sodium Assay can be used to discover and characterize the effects of many tens-of-thousands of compounds and environmental factors on effectors of intracellular Na<sup>+</sup>. ION's Brilliant Sodium Assay provides all the reagents necessary for use as a wash or no-wash assay with adherent or non-adherent cells. The optional use of a probenecid solution and an extracellular background masking solution (TRS) offers the ultimate in compatibility for cells types which are difficult to load with fluorescent Na<sup>+</sup> indicators (e.g. Chinese Hamster Ovary, CHO cells) and when performing assays in complete, serum-containing cell culture medium is desired.

ION's Brilliant Sodium Assay is compatible with fluorescence microscopes, flow cytometers, and plate readers capable of detecting fluorescein or more optimally, yellow fluorescent protein (YFP).

## Laboratory Procedures

### Getting Started

Before you begin, make sure that you have all the additional reagents and materials you will need for the successful completion of your experiment. While the ION Brilliant Sodium Assay Kit contains all the reagents you will need to prepare your cells for testing, your experiments will likely include other reagents which are not included in your ION Brilliant Sodium Assay Kit. Notably compounds to be tested are not included, neither are buffers and solvents for the dissolution of those compounds. The Brilliant Sodium Assay Kit also does not contain reagents necessary for cell culture. ION's pIONeer Kits do not contain DMSO for solubilizing the dye. In addition to reagents, a fluorescence plate reader that is capable of providing excitation between 485 - 525 nm and collecting emission at ~545 nm is required. Although kinetic plate readers capable of providing readouts at ~1 Hz, such as WaveFront Panoptic, Hamamatsu FDSS, Molecular Devices FLIPR and Molecular Devices FlexStation, are commonly used - ION's Brilliant Sodium Assay can also be used in an endpoint format on standard plate readers.

### Wash Method – Adherent Cells Only

The instructions given below are for one, 384-well microplate. Certain aspects of the instructions may need to be altered, as appropriate, for multiple microplates or other assay formats (e.g. 96-well microplates or non-adherent cells). The ING-2 AM and ING-2 AM-containing solutions should be protected from direct light.

1. Add 20 µL DMSO<sup>1</sup> (Reagent B) to the tube containing ING-2 AM (Reagent A).
2. Vortex until Reagent A is fully dissolved.
3. Add appropriate volume of water (Table 2) to a 15 mL centrifuge tube.
4. Add 1 mL of 10X Brilliant Assay Buffer (Reagent D) to tube from step 3.
5. Add 200 µL of DySolv (Reagent C) to the tube from step 4.
6. If desired add 200 µL of Probenecid Solution (Reagent F) to the tube from step 5.
7. Add 20 µL of ING-2 AM Solution from step 2 to the tube from step 6.
8. Briefly vortex the Dye Loading Solution, tube from step 7, to mix.

*Procedure Continues on Next Page*

Table 2		Dye Loading Solution (Wash Method)	
Label	Name	Method A	Method B
Reagent A+B	ING-2 AM Solution	20 µL	20 µL
Reagent C	50X DySolv	200 µL	200 µL
Reagent D	10X Brilliant Assay Buffer	1 mL	1 mL
Reagent F	50X Probenecid <sup>2</sup>	-	200 µL
	Water	8.6 mL	8.4 mL
	Total	10 mL	10 mL

9. Remove the cell-culture medium from the 384-well microplate containing the cells of interest.
10. Add 20  $\mu$ L per well of the Dye Loading Solution from step 8 to the microplate from step 9.
11. Incubate the microplate containing the cells and Dye Loading Solution for 30 minutes - 1 hour at 37° C.

**Table 3**      **Wash Solution**

Label	Name	Method A	Method B	Method C	Method D
Reagent D	10X Brilliant Assay Buffer	1 mL	1 mL	1 mL	1 mL
Reagent E	50X TRS <sup>2</sup>	-	200 $\mu$ L	-	200 $\mu$ L
Reagent F	50X Probenecid <sup>2</sup>	-	-	200 $\mu$ L	200 $\mu$ L
	Water	9 mL	8.8 mL	8.8 mL	8.6 mL
	Total	10 mL	10 mL	10 mL	10 mL

12. Prepare Wash Solution in a 15 mL centrifuge tube by adding the appropriate amounts of water, 10X Brilliant Assay Buffer (Reagent D) and other components if desired as shown in Table 3.
13. Briefly vortex the tube from step 12 to mix.
14. Remove Dye Loading Solution from microplate in step 11.
15. Add 20  $\mu$ L per well of the Wash Solution prepared in step 13 to the microplate from step 14.
16. Transfer the dye-loaded, cell-containing microplate from step 11, along with an additional microplate containing a stimulus solution of interest, to a kinetic-imaging plate reader (e.g. WaveFront Panoptic, Hamamatsu FDSS, Molecular Devices FLIPR or Molecular Devices FlexStation).
17. Acquire data using an excitation wavelength of  $\sim 520$  nm<sup>5</sup>, an emission wavelength of  $\sim 545$  nm and an acquisition frequency of  $\sim 1$  Hz.<sup>4</sup> Begin data acquisition and after 20 seconds add 5  $\mu$ L of the 5X stimulus solution to the cell-containing plate and continue data acquisition for an additional 90 seconds<sup>3</sup>.

## No-wash Method – Adherent Cells

1. Add 20  $\mu$ L DMSO<sup>1</sup> (Reagent B) to the tube containing ION ING-2 AM (Reagent A)
2. Vortex until Reagent A is fully dissolved.
3. Add appropriate volume of water (Table 4) to a 15 mL centrifuge tube.
4. Add 1 mL of 10X Brilliant Assay Buffer (Reagent D) to tube from step 3.
5. Add 200  $\mu$ L of DySolv (Reagent C) to the tube from step 4.
6. Add 200  $\mu$ L of TRS (Reagent E) to the tube from step 5.
7. If desired, add 200  $\mu$ L of Probenecid Solution (Reagent F) to the tube from step 6.
8. Add 20  $\mu$ L of ING-2 AM Solution from step 2 to the tube from step 7.
9. Briefly vortex the Dye Loading Solution, tube from step 8, to mix.

*Procedure Continues on Next Page*

**Table 4** Dye Loading Solution (No-wash Method)

Label	Name	Method A	Method B
Reagent A+B	ING-2 AM Solution	20 µL	20 µL
Reagent C	50X DySolv	200 µL	200 µL
Reagent D	10X Brilliant Assay Buffer	1 mL	1 mL
Reagent E	50X TRS <sup>2</sup>	200 µL	200 µL
Reagent F	50X Probenecid <sup>2</sup>	-	200 µL
	Water	8.6 mL	8.4 mL
	Total	10 mL	10 mL

10. 13. Remove the cell-culture medium from the 384-well microplate containing the cells of interest.
11. Add 20 µL per well of the Dye Loading Solution from step 9 to the cell-containing microplate.
12. Incubate the microplate containing the cells and Dye Loading Solution for 30 minutes - 1 hour at 37°C in a cell culture incubator.
13. Transfer the dye-loaded, cell-containing microplate from step 11, along with an additional microplate containing a stimulus solution of interest, to a kinetic-imaging plate reader (e.g. WaveFront Panoptic, Hamamatsu FDSS, Molecular Devices FLIPR or Molecular Devices FlexStation). See Table 5 for recommended instrument settings.
14. Acquire data using an excitation wavelength of ~ 515 nm<sup>5</sup>, an emission wavelength of ~ 545 nm and an acquisition frequency of ~1 Hz.<sup>4</sup> Begin data acquisition and after 20 seconds add 10 µL of the 5X stimulus solution to the cell-containing plate and continue data acquisition for an additional 90 seconds<sup>3</sup>.

**Table 5** Recommended Instrument Settings

Setting	Recommendation
Read Mode	'Bottom' read mode only
Ex/Em wavelengths <sup>4</sup>	~515 nm/545 nm
Cutoff wavelength	530 nm
Filter selection	GFP, FITC, or YFP

Contact [support@ionbiosciences.com](mailto:support@ionbiosciences.com) for additional recommendations and guidance on optimizing your application.

## Laboratory Procedures - Footnotes

- <sup>1</sup> DMSO is hygroscopic and should be stored tightly closed at -20° C with desiccant pack to prevent the solvent from becoming wet. Wet solvent causes difficulties with dissolution of the dye. Use the DMSO within 6 months of receipt. ION's pIONeer Kits do not contain DMSO.
- <sup>2</sup> Caution is advised when using Probenecid and/or TRS as they may have undesirable effects on assay performance for the target of interest. Probenecid may be included in the Dye Loading Solution to aid in dye retention. This may be particularly important in certain cell lines (e.g. CHO cells). TRS contains a membrane-impermeant dye useful for masking extracellular fluorescence.
- <sup>3</sup> The timing and volume of stimulus solution addition may vary. Some experiments may include the addition of other solutions to the cell-containing microplate prior to the addition of the stimulus solution. In these cases, the volume of the stimulus solution addition should be altered to account for the additional volume of solution in the cell-containing microplate.
- <sup>4</sup> For targets where changes in intracellular sodium concentrations are slow or sustained, an endpoint assay format can be used. We recommend acquiring data before the addition of stimulus ( $F_0$ ) and again 15-30 min after the addition of stimulus.
- <sup>5</sup> Excitation sources commonly used for fluorescein, Fluo-4, and GFP (480–490 nm) are also compatible with Brilliant Sodium.
- <sup>6</sup> To prevent bleed-through or spectral overlap, the Ex/Em wavelengths may need to be optimized by broadening the interval between the wavelengths. Results generated on a Molecular Devices Flexstation 3.

## Additional Information

Additional dye indicator and buffer reagents can be purchased either directly from our website or by contacting our Sales Department. Custom and bulk sizes are also available. Contact Sales for more information.

Table 6		Additional Reagents		Available Sizes	
Kit Label	Name	Size		Catalog #	
Reagent A	ING-2 AM	500 µg x 1 Vial		2011C	
		50 µg x 10 Vials		2011F	
		50 µg x 3 Vials		2011G	
Reagent C	50X DySolv	20 mL Bottle		7501A	
Reagent D	10X Brilliant Assay Buffer	10 mL Bottle		7010X	
Reagent E	50X TRS	20 mL Bottle		7060A	
Reagent F	50X Probenecid	20 mL Bottle		7300P-50	