

HEK293-huK_{ir}3.1/3.2 Recombinant cell line

Catalog No. C1101

Introduction

G protein-gated inwardly rectifying potassium channels (GIRK) are a large family of inwardly rectifying potassium channels that play important roles in neurophysiology and cardiac physiology. Mammals express four GIRK subunits (K_{ir}3.1-K_{ir}3.4); the predominant form of GIRK in the brain is a K_{ir}3.1/K_{ir}3.2 heterotetramer. GIRK channels are opened by the direct binding of the G protein $\beta\gamma$ subunit (G $\beta\gamma$), which is released from the heterotrimeric G protein (G $\alpha\beta\gamma$) upon the activation of G protein coupled receptors (GPCRs). GIRK is exclusively associated with G_{i/o} GPCRs such as dopamine receptors, serotonin receptors, muscarinic receptors and opioid receptors. Due to its critical role of regulating cell excitability, GIRK presents a potential therapeutic target for epilepsy, Parkinson's disease, Down's syndrome, pain, drug addiction, and many more pathologies.

Recombinant HEK293-huK_{ir}3.1/3.2 cell line express both human KCNJ3 (K_{ir}3.1, potassium inwardly rectifying channel subfamily J member 3, accession number NM_002239) and human KCNJ6 (K_{ir}3.2, potassium inwardly rectifying channel subfamily J member 6, accession number NM_002240). K_{ir}3.1 is unable to independently form functional channels on the plasma membrane; it forms heterotetrameric channels. By contrast, K_{ir}3.2 can form homotetramers and can assemble with other GIRK subunits to form heterotetramers, such as K_{ir}3.1/K_{ir}3.2. K_{ir}3.1/K_{ir}3.2 channels have a greater tendency to allow potassium flow into a cell than out of a cell; they modulate many physiological processes, including heart rate in cardiac cells and circuit activity in neuronal cells. Our HEK293-huK_{ir}3.1/3.2 cell line is suitable for discovering modulators of GIRK channel activity and screening for agonist/antagonist of various G_{i/o} GPCRs that associate with GIRK channels.

Materials Provided

One vial of HEK293-huK_{ir}3.1/3.2 recombinant cells, 2 x 10⁶ cells in 1 mL of Bambanker serum free freezing medium.

Storage

Cells are shipped on dry ice and should arrive frozen. To ensure maximum cell viability, store the cell vial in liquid nitrogen immediately upon receipt.

Mycoplasma testing

The cell line has been screened using ATCC PCR-based testing service which covers 60 species of *Mycoplasma*, *Acholeplasma*, *Spiroplasma* and *Ureaplasma* including the eight species most likely to afflict cell cultures: *M. arginini*, *M. fermentans*, *M. hominis*, *M. hyorhina*, *M. orale*, *M. pirum*, *M. salivarium*, and *A. laidlawii*. The absence of Mycoplasma species has been confirmed for each lot. For detailed verification, please see the product certificate of analysis.

Materials Required But Not Supplied

| Category | Name | Recommendations |
|--|---------------------------|---|
| Reagents | Cell culture media | Advanced MEM (ThermoFisher Cat#: 12491013) |
| | GlutaMax™ | GlutaMAX™ Supplement (ThermoFisher Cat#: 35050061) |
| | Fetal bovine serum (FBS) | Fetal Bovine Serum (Corning Cat#: 35-016-CV) |
| | Blasticidin | Blasticidin S HCl (ThermoFisher Cat#: A1113903) |
| | Puromycin | Puromycin Dihydrochloride (ThermoFisher Cat#: A1113803) |
| | TrypLE™ | TrypLE™ Express Enzyme (ThermoFisher Cat#: 12605028) |
| | Phosphate buffered saline | PBS, pH 7.4 (ThermoFisher Cat#: 10010023) |
| | Freezing media | Bambanker Serum Free Cell Freezing Medium |
| | Equipment | Single and multichannel micropipettes and pipette tips |
| 50 mL and 15 mL conical centrifuge tubes | | |
| 1.5 mL microtubes | | |
| Tissue culture disposable pipettes (1 mL-25 mL) and tissue culture flasks (e.g. T75 flask) | | |
| 0.2 µm filter unit for medium sterilization | | |
| Cryovials for freezing cells | | |
| Automated Cell Counter (or Hemocytometer) | | |
| Humidified tissue culture incubator (37°C and 5% CO ₂) | | |

Cell Culture Protocol

A. Cell culture media

1. Make thawing medium by combining 25 mL of FBS, 5 mL of 100X GlutaMax with 470 mL 1X Advance MEM media for the final concentration of 5% FBS and 1X GlutaMax.
2. To make complete culture medium, add 250 µL of 10 mg/mL Blasticidin, and 150 µL of 10 mg/mL Puromycin into 500 mL of medium from **step 1** for the final concentrations of 5 µg/mL Blasticidin and 3 µg/mL Puromycin.
3. (Optional) Sterilize all media using 0.2 µm filter.
4. Pre-warm all media in a 37°C water bath prior to use.

B. Thawing cells

1. Aliquot 8 mL of thawing medium in a 15 mL conical centrifuge tube.
2. Place the cryovial in a 37°C water bath briefly, until only small ice crystals remain and the cell pellet is almost completely thawed. The thawing time typically ranges from 2 to 3 minutes. DO NOT vortex freshly thawed cells.
3. Spray and wipe the external surface of the cryovial with 70% ethanol. Transfer the vial to a biosafety cabinet.

Cell Culture Protocol Cont.

4. Gently add the cells into the pre-filled 15 mL conical tube dropwise.
5. Rinse the cryovial with 1 mL of cell suspension to maximize cell recovery and add it back into the 15 mL conical tube.
6. Centrifuge the 15 mL conical tube containing cells at 500xg for 5 min at 25°C.
7. Add 14 mL of complete cell culture medium into a T75 flask.
8. After centrifugation, carefully aspirate the supernatant without disturbing the cell pellet.
9. Gently resuspend the cell pellet with 1 mL of complete cell culture medium.
10. Transfer the cell suspension to the T75 flask and gently swirl the flask to distribute the cells evenly in solution.
11. Incubate the flask in a 37°C and 5% CO₂ humidified cell culture incubator.
12. Maintain the cells in culture by exchanging culture media every 2-3 days, until they reach >80% confluency in a T75.

C. Cell passage

1. Remove the T75 flask from the tissue culture incubator and place it in a biosafety cabinet.
2. Gently aspirate the media from the T75 flask.
3. Add 15 mL of PBS into the T75 flask, and gently rock the flask back and forth to rinse the cells.
4. Gently aspirate PBS from the flask.
5. Add 3 mL of pre-warmed TrypLE to the flask. Gently rock the flask back and forth to ensure that the flask is uniformly covered with TrypLE.
6. Incubate the flask at 37°C and 5% CO₂ until the cells have detached. Usually it will take 8 - 10 min.
7. Remove the flask from the incubator and observe the cells under a microscope to confirm cell detachment. If necessary, gently tap the edge of the flask to detach cells from the surface.
8. Add 10 mL of pre-warmed complete cell culture media to the detached cells in the flask.
9. Transfer cell suspension into a 15 mL conical tube.
10. Centrifuge the 15 mL conical tube containing cells at 500xg for 5 min at 25°C.
11. After centrifugation, carefully aspirate the supernatant without disturb the cell pellet.
12. Resuspend the cell pellet with 5 mL of complete cell culture medium. Pipette the cell suspension up and down several times to generate a single cell suspension without any clumps.
13. Take a small aliquot of the cell suspension (~20 µL) for cell counting, and determine the number of cells and volume required to seed into a T75 flask. HEK293-huK_i3.1/3.2 cell doubling time is ~36h. A confluent T75 flask normally yields ~10x10⁶ cells. We recommend to seed 1 X 10⁶ cells for a T75 flask.
14. Add 15 mL of complete cell culture medium into a new T75 flask, followed by addition of the appropriate volume of cell suspension. Transfer the flask to a tissue culture incubator, and incubate the cells at 37°C and 5%CO₂.

Cell Culture Protocol Cont.

D. Cryopreservation

1. Harvest cells according to the description in **section C Cell Passage step 1 - 9**.
2. Set aside a small fraction of the suspended cells (100 μ L or less) in a separate tube for cell counting.
3. Count cells, calculate the concentration of cells and the total number of cells in the original 15 mL conical tube.
4. Centrifuge the 15 mL conical tube containing cells at 500xg for 5 min at 25°C.
5. After centrifugation, discard the supernatant, being careful not to disturb the cell pellet.
6. Based on the total cell number calculated in **Step 3**, resuspend the cells to the desired concentration (e.g. 2.0×10^6 cells/mL) in Bamberker serum free cell freezing medium.
7. Aliquot 1 mL of the cell suspension into each of the labeled 2 mL cryovials. Seal the cryovials tightly.
8. Freeze cells in a -80°C freezer at a controlled rate of -1°C/minute overnight in a commercially-available freezing container.
9. The following day, transfer the vials into the vapor phase of liquid nitrogen for long-term storage.

Functional Validation

The function of $K_{ir}3.1/3.2$ channels were validated by GIRK potassium channel activator ML297 (VU0456810) (Figure 1) and VU551 (VU0529331) (Figure 2). In addition, $G_{i/o}$ GPCRs were co-expressed into HEK293-hu $K_{ir}3.1/3.2$ cells; GPCR agonists activated $K_{ir}3.1/3.2$ channels, resulting in an increased Tl^+ influx detected using $G_{i/o}$ GPCR-GIRK thallium flux assay (Figure 3). HEK293-hu $K_{ir}3.1/3.2$ cells co-expressing dopamine D2 receptors showed distinct signal increases upon activation by dopamine in a HTS assay with Z' factor = 0.55 (Figure 4).




Figure 1. ML297 Dose Response Curve. Thallium flux assay was performed using FlexStation 3 kinetic plate reader with excitation wavelength 490 nm and emission wavelength 520 nm. HEK293- $K_{ir}3.1/3.2$ cells were treated with potassium channel activator ML297 at the highest concentration of 10 μ M. Dose response curves generated an EC_{50} : 100 nM.




Figure 2. VU551 Dose Response Curve. Thallium flux assay was performed using FlexStation 3 kinetic plate reader with excitation wavelength 490 nm and emission wavelength 520 nm. HEK293- $K_{ir}3.1/3.2$ cells were treated with potassium channel activator VU551 at the highest concentration of 10 μ M. Dose response curves generated an EC_{50} : 143 nM.

Functional Validation Count.




Figure 3. 96-well $G_{i/o}$ GPCR-GIRK Thallium Flux Assay Dose Response Curves. Dopamine D2 receptor, Serotonin 1A receptor, Muscarinic M2 receptor, and delta opioid receptor were transduced in HEK293-huK_i3.1/3.2 cells using the optimized amount of transduction conditions for 24h. Thallium flux assay was performed using a FlexStation 3 kinetic plate reader with excitation wavelength 490 nm and emission wavelength 520 nm. Dopamine, Serotonin, Acetylcholine, and Enkephalin were tested with the highest concentration of 10 μ M. Dose response curves generated dopamine EC₅₀: 6.7 nM, serotonin EC₅₀: 47.8 nM, acetylcholine EC₅₀: 15.3 nM, and Enkephalin EC₅₀: 0.84 nM.




Figure 4. 96-well $G_{i/o}$ GPCR-GIRK Thallium Flux Assay Z' Factor Analysis. Dopamine D2 receptor was transduced in HEK293-huK_i3.1/3.2 cells using the optimized transduction conditions for 24h. Thallium Flux assay was performed using FlexStation 3 kinetic plate reader with excitation wavelength 490 nm and emission wavelength 520 nm. For positive control group 10 μ M of dopamine was used to activate dopamine D2 receptor; for negative control group buffer alone was added to the plate. Normalized kinetic data showed a clear distinction between positive and negative control groups. The calculated Z' factor is 0.55 indicating a high quality screening assay.

Vector and Sequence

Vector description

1. Mammalian Gene Expression Vector, pRP[Exp]-Bsd-CAG>hKCNJ3[NM_002239.4]
2. Mammalian Gene Expression Vector, pPB[Exp]-Puro-CAG>hKCNJ6[NM_002240.5]

hKCNJ3 sequence (NM_022039.4)

MSALRRKFGDDYQVVTSSSGSLQPQPGQDPQQQLVPKKRQRFVDKNGRCNVQHGNLGSETSRYLSDLFTTLVLDLKWVWNLFI
TYTVAWLFMASMWWVIAAYTRGDLNKAHVGNYPVCVANVYFNSAFLFFIETEATIGYGYRYITDKCPEGIILFLFQSILGSIVDAFLIGCMFI
KMSQPKKRAETLMFSEHAVISMVDGKLTLMFRVGNLRNSHMVSAQIRCKLLKSRQTPEGEFLPLDQLELDVGFSTGADQLFLVSPLTICH
VIDAKSPFYDLSQRSMQTEQFEIVVILEGIVETTGMTQARTSYTEDEVLWGHRRFPVISLEEGFFKVDYSQFHATFEVPTPPYSVKEQEE
MLLMSSPLIAPAITNSKERHNSVECLDGLDDITTKLPSKLQKITGREDFPKLLRMSSTTSEKAYSLGDLPMKLQRISSVPGNSEKLVSKT
TKMLSDPMSQSVADLPPKLQKMAGGAARMEGNLPAKLRKMNSDRFT

hKCNJ6 sequence (NM_002240.5)

MAKLTESMTNVLEGDSMDQDVESPAIHQPPLPKQARDDLPRHISRDRTRKRIQRYVRKDGKCNVHHGNVRETYRYLTDIFTTLVLDLKW
RFNLLIFVMVYTVTWLFFGMWWLIAYIRGDMHDIEDPSWTPCVTNLNGFVSAFLFSIETETTIGYGYRVITDKCPEGIILLIQSVLGSIVNA
FMVGCMEFVKISQPKKRAETLVFSTHAVISMVDGKLTLMFRVGNLRNSHIVEASIRAKLIKSQTSEGEFIPLNQTDINVGYYTGDDRLFLVS
PLIISHEINQQSPFWEISKAQLPKEELEIVVILEGMVEATGMTQARSSYITSEILWGYRFTPVLTEDEGFVEVDYNSFHETYETSTPPLSAK
ELAELASRAELPLSWSVSSKLNQHAETEETEEENLEEQTERNGDVANLENESKV

References

1. Luscher C, Slesinger PA. [Emerging roles for G protein-gated inwardly rectifying potassium \(GIRK\) channels in health and disease](#). *Nat. Rev. Neurosci.* **2010**, 11(5):301-315.
2. Kano HK, Toyama Y, Imai S, Iwagashi Y, Mase Y, Yokogawa M, Osawa M, Shimada I. [Structure mechanism underlying G protein family-specific regulation of G protein-gated inwardly rectifying potassium channel](#). *Nat. Commun.* **2019**, 10(1):2008.
3. Liao YJ, Jan YN, Jan LY. [Heteromultimerization of G-protein-gated inwardly rectifying K⁺ channel proteins GIRK1 and GIRK2 and their altered expression in weaver brain](#). *J. Neurosci.* **1996**, 16, 7137-7150.
4. Kaufmann K, Romaine I, Days E, Pascual C, Malik A, Yang L, Zou B, Du Y, Sliwoski G, Morrison RD, Denton J, Niswender CM, Daniels JS, Sulikowski GA, Xie XS, Lindsley CW, Weaver CD. [ML297 \(VU0456810\), the First Potent and Selective Activator of the GIRK Potassium Channel, Displays Antiepileptic Properties in Mice](#). *ACS Chem. Neurosci.* **2013**, 4(9): 1278-1286.
5. Kozek KA, Du Y, Sharma S, Prael FJ, Spitznagel BD, Kharade SV, Denton JS, Hopkins CR, Weaver CD. [Discovery and Characterization of VU0529331, a Synthetic Small-Molecule Activator of Homomeric G Protein-gated, Inwardly-rectifying, Potassium \(GIRK\) Channels](#). *ACS Chem. Neurosci.* **2019**, 10(1): 358-370.
6. Tateyama M, Kubo Y. [Gi/o-coupled muscarinic receptors co-localized with GIRK channel for efficient channel activation](#). *PLoS One.* **2018**, 13(9):e0204447.
7. Gadgaard C, Jensen AA. [Functional characterization of 5-HT_{1A} and 5HT_{1B} serotonin receptor signaling through G-protein-activated inwardly rectifying K⁺ channels in a fluorescence-based membrane potential assay](#). *Biochem. Pharmacol.* **2020**, 175:113870.

Limited Use License Agreement

Cell Line Terms and Conditions are explicitly accepted upon purchase. An additional agreement is needed for commercial use. These cells are sold under a non-transferable, limited use license to use the cells and their progeny in your immediate laboratory, for research use only. Commercial use of the cells is not permitted under this license. This license does not permit you to share, transfer, distribute, sell, sublicense, or otherwise make the cells available for use to other researchers, laboratories, departments, research institutions, hospitals, universities, CROs, or biotech companies. The license does not permit use of the cells in humans or for therapeutic or drug use. Inappropriate use or distribution of the cells will result in revocation of this license and result in an immediate cease of sales and distribution of ION Biosciences products to your laboratory.

ION Biosciences does not warrant the suitability of the cells for any particular use, and does not accept any liability including direct, indirect, consequential, or incidental damages in connection with the handling or use of the cells. Modifications of the cells not discussed herein, transfer to another facility, or commercial use requires a separate license and additional fees. Commercial use means any activity by a party for consideration and may include, but is not limited to (1) use of the product or associated components to provide a service, information or data to a third party and/or (2) resale of the product or associated components, whether or not such product or its components are resold for use in research. Contact sales@ionbiosciences.com for details.