

Thallium flux assay adaptation for multi-instrument compatibility

Derek S. Hernandez, C. David Weaver
ION Biosciences, San Marcos, TX USA



Abstract

Over the past 15+ years, fluorescence-based measures of Tl^+ flux have brought about the discovery of small-molecule modulators of a host of ion channels, transporters, GPCRs and other targets of interest for both drug discovery and basic research.^{1,2} Here, we introduce ION's **Brilliant Thallium Snapshot** Assay Kit, which provides a brand new way to conduct thallium flux assays. **Brilliant Thallium Snapshot** is designed for multi-well plate-based, high-throughput measurements of Tl^+ flux through K^+ , Na^+ , non-selective cation channels, and some Na^+ or K^+ transporters. Unlike traditional thallium flux assays, our **patent-pending Snapshot** assay format generates a **long-lasting signal** so you can easily detect and quantify changes in ion channel or transporter activity on most fluorescence-capable instruments. So whether you have a FLIPR® (HTS), a fluorescence microscope, a standard fluorescence plate reader, a high-content imager (HCS), or a flow cytometer - you can run functional screens on viable cells using thallium flux, enabling entirely new possibilities.

Instrument compatibility is demonstrated using a CHO GIRK1/2 (CHO G12) expressing cell line. GIRK channels are modulated using known activators, ML297³ or VU0466551 (VU551).⁴ Comparable EC_{50} values are obtained using data acquired from three commonly available instruments -



a fluorescence microscope, standard fluorescence plate reader, and Flexstation®. Thallium flux paired with microscopy enables monitoring of specific cells within a diverse population - an ideal solution for measuring cell-specific ion channel activity in complex co-culture systems. Other applications of this technology could include clonal selection using FACS and identifying target expressing cells within dissociated tissue samples.



Methods

A simple workflow designed to accommodate a variety of instruments, automation, and users at all skill levels.

- 1) Prepare cell containing microplate
- 2) Remove medium (if possible)
- 3) Load cells with **dye loading solution** (contains Thallos AM and Tl^+), and incubate for 30 - 60 minutes
- 4) Acquire initial fluorescence readings (F_0) or images using appropriate filters (Ex/Em: 485/525 nm)
- 5) Add compound(s) of interest, and incubate for 5 - 30 minutes
- 6) Acquire fluorescence readings or images



Results

Figure 1. Traditional thallium flux assay. Rapid Tl^+ influx kinetics and fast equilibration limit the signal detection window. Inset) VU551 dose response curve with reported $EC_{50} = 77$ nM. Denoted concentrations are in μM .





Figure 2. Kinetic fluorescence profiles of CHO GIRK1/2 (G12) cells using Brilliant Thallium Snapshot. A) Cells are pre-loaded with Thallos AM and Tl^+ prior to adding the GIRK activator, VU551. Kinetic data acquired using a Flexstation®. Efflux of Tl^+ upon GIRK channel activation yields a decrease in fluorescence. B) Changes in signal are sustained for >40 minutes as demonstrated using a standard fluorescence plate reader (Cytation 5). Reads were acquired at 1 min intervals. Denoted VU551 concentrations are in μM .




Figure 3. Brilliant Thallium Snapshot assays on a plate reader and microscope. A) Checkerboard layout of a 96-well plate was used to calculate a Z' . Vehicle and activator (1 μM VU551) treated wells are easily distinguished 20 minutes after compound addition ($Z' = 0.56$). B) Representative fluorescence images of vehicle- (left) and VU551-treated (right) CHO G12 cells demonstrating a 43% reduction in mean cellular fluorescence. Images acquired 20 min after VU551 addition using fluorescein filters and a 4X objective.





Figure 4. Dose-response curves of CHO G12 cells to VU551 measured using several instruments. When using Brilliant Thallium Snapshot Assay, comparable EC_{50} values are obtained on various fluorescence-capable instruments. All values are similar to data acquired using ION's Brilliant Thallium Flex kit ($EC_{50} = 77$ nM) and published literature ($EC_{50} = 75$ nM).

Results




Figure 5. Imaging-based GIRK channel activity assays in a co-culture system. A) CHO K1 (pink) cells were stained with CytoTracker Red prior to plating in the same wells as CHO cells overexpressing GIRK1/2 (CHO G12 - yellow). CytoTracker Red fluorescence was used to identify cell type. B) Overlay of GFP and TRITC images acquired before the addition of VU551 (F_0). C) Overlay of GFP and TRITC images acquired 10 min after the addition of VU551 (F_{10}). A significant decrease in fluorescence was observed only in CHO G12 cells in the presence of VU551. Images (A - C) display the same field of view.




Figure 6. Flow cytometry-based GIRK channel assays. Left) Median fluorescence of separate cell populations acquired before and at reported timepoints after the addition of HHBSS or 3 μM VU551. Right) Histograms acquired 50 minutes post-treatment reveal a significant drop in cellular fluorescence in the group treated with a GIRK activator (1 μM VU551) relative to a control (HHBSS).

Conclusions

ION's **Brilliant Thallium Snapshot** assay delivers:

- Sustained signal amplitude
- Increased accessibility
- Cell-specific readouts
- Greater convenience
- Reduced interference from background flux

For additional information, visit product page here:



References

1. Weaver CD. Thallium Flux Assay for Measuring the Activity of Monovalent Cation Channels and Transporters. *Methods Mol Biol.* 2018.
2. Niswender CM, Johnson KA, Luo Q, et al. A novel assay of Gi/o-linked G protein-coupled receptor coupling to potassium channels provides new insights into the pharmacology of the group III metabotropic glutamate receptors. *Mol Pharmacol.* 2008.
3. Kaufmann K, Romaine I, Days E, et al. ML297 (VU0456810), the first potent and selective activator of the GIRK potassium channel, displays antiepileptic properties in mice. *ACS Chem Neurosci.* 2013.
4. Wen W, Wu W, Weaver CD, Lindsley CW. Discovery of potent and selective GIRK1/2 modulators via 'molecular switches' within a series of 1-(3-cyclopropyl-1-phenyl-1H-pyrazol-5-yl)ureas. *Bioorg Med Chem Lett.* 2014.