

THALLOS GOLD AM

Table 1 - Materials Needed		Most Items Available from ION Biosciences	
Name	Volume	Description/Purpose	
Thallos Gold AM (50 µg vial)	Dry - 1 Vial	Membrane Permeable Thallium Dye Indicator	
DMSO ¹	25 µL	Solvent for Dissolution of Dye	
100X Pluronic F-127 solution ²	100 µL	Biocompatible Surfactant for Dye Loading	
100X Probenecid solution (optional) ³	100 µL	Intracellular Dye Retention Aid	
50X TRS (optional) ⁴	200 µL	Extracellular Fluorescence Masking Agent	
Assay Buffer - We Suggest One of the Following:			
⇒ 1X HEPES-Buffered Hank's Balanced Salt Solution	10 mL	1X - Ready to Use Assay Buffer	
⇒ 10X Brilliant Thallium Assay Buffer	1 mL	10X - Concentrated Assay Buffer	
10X Brilliant Chloride-Free Stimulus Buffer	0.5 - 1 mL	Stimulus Component [Na ⁺]=1.38 M [K ⁺]=0 M	
10X High-Potassium, Chloride-Free Stimulus Buffer	0.5 mL	Stimulus Component [Na ⁺]=0 M [K ⁺]=1.38 M	
50 mM Thallium Sulfate (Tl ₂ SO ₄) Solution	0.5 mL	Stimulus Component - Thallium Ion (Tl ⁺)	
Water	As Needed	Dilution of 10X Buffer (if used) and Stimulus	

Description

Thallos Gold is a gold fluorescent (Ex/Em = 528/552 nm) intracellular thallium (Tl⁺) indicator. Thallos Gold delivers outstanding results for fluorescence-based potassium (K⁺) channel HTS and for a wide variety of monovalent cation (sodium - Na⁺) channels, transporters, and GPCRs. Red-shifted excitation and emission enable multiplexing with GFP-expressing cells or other green fluorescent indicators and minimizes interference from auto-fluorescent compounds.

Laboratory Procedures

The following protocol provides general guidelines for using this dye to measure intracellular thallium flux. All loading conditions (dye concentration, temperature, and time) should be optimized for your specific assay, application, and instrumentation.

1. Allow all reagents to warm to room temperature before proceeding.
2. Add 25 µL DMSO¹ to the tube containing Thallos Gold AM.
3. Vortex until Thallos Gold AM is fully dissolved. Centrifuge briefly to collect all contents at the tube bottom.

Procedure Continues on Next Page

Laboratory Procedures (continued)

4. Add the appropriate volume (see **Table 2**, below) of 1X Assay Buffer to a 15 mL conical tube. If using 10X Assay Buffer, add the appropriate volume of Water (**Table 3**, below) and 1 mL of 10X Assay Buffer to a 15 mL conical tube.
5. Add 100 µL of 100X Pluronic F-127² solution to the conical tube from **step 4**.
6. (Optional) Add 100 µL of 100X Probenecid³ solution to the conical tube from **step 5**.
7. (Optional) Add 200 µL of 50X TRS⁴ solution to the conical tube from **step 6**.
8. Vortex conical tube from **step 7** briefly to mix.
9. Add the entire contents of the Thallos Gold AM in DMSO solution from **step 3** to the conical tube from **step 8** to make the **Dye Loading Solution**⁵.
10. Vortex the **Dye Loading Solution**⁵ from **step 9** briefly to mix.

Table 2	Dye Loading Solution ⁵				Using 1X Assay Buffer			
Name	Method A	Method B	Method C	Method D	Method A	Method B	Method C	Method D
Thallos Gold AM in DMSO ¹ Solution	25 µL	25 µL	25 µL	25 µL	25 µL	25 µL	25 µL	25 µL
100X Pluronic F-127 ² solution	100 µL	100 µL	100 µL	100 µL	100 µL	100 µL	100 µL	100 µL
100X Probenecid ³ solution (optional)	100 µL	100 µL	-	-	100 µL	100 µL	-	-
50X TRS ⁴ (optional)	200 µL	-	200 µL	-	200 µL	-	200 µL	-
1X HEPES-Buffered Hanks Balanced Salt Solution	9.6 mL	9.8 mL	9.7 mL	9.9 mL	9.6 mL	9.8 mL	9.7 mL	9.9 mL
Total	10 mL	10 mL	10 mL	10 mL	10 mL	10 mL	10 mL	10 mL

Table 3	Dye Loading Solution ⁵				Using 10X Assay Buffer			
Name	Method A	Method B	Method C	Method D	Method A	Method B	Method C	Method D
Thallos Gold AM in DMSO ¹ Solution	25 µL	25 µL	25 µL	25 µL	25 µL	25 µL	25 µL	25 µL
100X Pluronic F-127 ² solution	100 µL	100 µL	100 µL	100 µL	100 µL	100 µL	100 µL	100 µL
100X Probenecid ³ solution (optional)	100 µL	100 µL	-	-	100 µL	100 µL	-	-
50X TRS ⁴ (optional)	200 µL	-	200 µL	-	200 µL	-	200 µL	-
10X Brilliant Thallium Assay Buffer	1 mL	1 mL	1 mL	1 mL	1 mL	1 mL	1 mL	1 mL
Water	8.6 mL	8.8 mL	8.7 mL	8.9 mL	8.6 mL	8.8 mL	8.7 mL	8.9 mL
Total	10 mL	10 mL	10 mL	10 mL	10 mL	10 mL	10 mL	10 mL

Procedure Continues on Next Page

Laboratory Procedures (continued)

11. Remove the cell culture medium and add the **Dye Loading Solution**⁵ from **step 10**. Recommend volumes are: 35 mm dish or 6-well plate, 1.5 mL; 96-well plate, 100 μ L; 384-well plate, 20 μ L.⁶
12. Incubate in a cell culture incubator at 37°C for 30 - 60 minutes.
13. Prepare a Thallium Stimulus Solution containing ~ 5 mM of thallium (Tl^+) ion. **Table 4** below provides two examples of Thallium Stimulus solutions useful for many types of non-voltage-gated and voltage-gated monovalent cation channels and transporters. The concentration of thallium in the stimulus solution may be varied to achieve the desired result. The final thallium concentration in the cell-containing microplate wells post-thallium stimulus buffer addition should not exceed 4.8 mM due to the solubility limit of thallium in chloride-containing solutions.

Table 4 Thallium Stimulus Solution

Name	Example A	Example B
10X Chloride-Free Stimulus Buffer ($[K^+] = 0$ M, $[Na^+] = 1.38$ M)	1 mL	0.5 mL
10X High-Potassium, Chloride-Free Stimulus Buffer ($[K^+] = 1.38$ M, $[Na^+] = 0$ M)	-	0.5 mL
50 mM Thallium Sulfate (Tl_2SO_4) Solution	0.5 mL	0.5 mL
Water	8.5 mL	8.5 mL
Total	10 mL	10 mL

14. Conduct a wash⁶ step to remove the **Dye Loading Solution**⁵ and replace with cell culture medium or assay buffer⁷. Repeat, if necessary, to completely remove extracellular dye.
15. Acquire data using a fluorescence microscope⁷ equipped with YFP filters or a fluorescence plate reader using an excitation wavelength of ~525 nm, an emission wavelength of ~560 nm and an acquisition frequency of 1-10 Hz⁸. See **Table 5** below for recommended settings. Begin data acquisition and after 10 seconds, add 5 μ L (for a 384-well plate) or 25 μ L (for a 96-well plate) of the Thallium Stimulus Solution to each well of the cell-containing plate and continue data acquisition for an additional 90 seconds

Table 5 Recommended Instrument Settings

Setting	Recommendation
Read Mode (Plate Readers)	'Bottom' read mode only
Ex/Em wavelengths ⁸	~525 nm/560 nm
Cutoff wavelength	550 nm
Filter selection	YFP

Contact support@ionbiosciences.com for additional recommendations and guidance on optimizing to your application.

Laboratory Procedures - Footnotes

- ¹ DMSO is hygroscopic and should be stored tightly closed. Wet solvent causes difficulties with dissolution of the dye.
- ² Pluronic F-127 is a biocompatible surfactant that aids in dye dissolution, ensuring equitable dye distribution and cellular loading.
- ³ Probenecid is an anion transport inhibitor that improves intracellular dye retention. Although it is not required for all cell types and dyes, it is recommended in most cases to optimize assay performance.
- ⁴ TRS is a membrane impermeant dye useful for masking extracellular fluorescence. Caution is advised when using TRS or any other extracellular masking solutions as they may have undesirable effects on assay performance for the target of interest.
- ⁵ The Dye Loading Solution should be used within 2 hours of dye addition for best results.
- ⁶ In some cases, a no wash format works best. If a no wash format is indicated for your application, we recommend doubling the concentration of 100X Pluronic F-127, 100X Probenecid, and 50X TRS in your dye loading buffer.
- ⁷ To minimize extracellular background, the dye loading solution can be replaced with assay buffer containing 1X Probenecid solution and/or 1X TRS solution.
- ⁸ To prevent bleed-through or spectral overlap, the Ex/Em wavelengths may need to be optimized by broadening the interval between the wavelengths.

Additional Information

Dye indicator and buffer reagents can be purchased either directly from our website or by contacting our Sales Department. Custom and bulk sizes are also available. Contact Sales for more information.

Table 4	Available Reagents	Available Sizes	
Name		Size	Catalog #
Thallos Gold AM		500 µg x 1 Vial	1391C
		50 µg x 10 Vials	1391F
		50 µg x 3 Vials	1391G
100X Pluronic F-127 Solution		10 mL Bottle	7601A
100X Probenecid Solution		10 mL Bottle	7300P-100
50X TRS Solution		20 mL Bottle	7060A
1X HEPES-Buffered Hank's Balanced Salt Solution (1X HHBSS)		100 mL Bottle	7001
10X Brilliant Thallium Assay Buffer		10 mL Bottle	7010T
10X Brilliant Chloride-Free Stimulus Buffer		10 mL Bottle	7020B
10X High-Potassium, Chloride Free Stimulus Buffer		10 mL Bottle	7030S
50 mM Thallium Sulfate Solution		5 mL Bottle	7040S