

100X Pluronic F-127

Catalog # 7601A

Category	Specification
Product Details	
Intended Use	A solubilization aid for water insoluble dyes, used as part of a dye-loading solution
Product Size	10 mL (Contact Sales for custom and bulk sized orders)
Product Concentration	100X
Storage Temperature	4°C
Product Form	
Appearance	Clear, Colorless solution
Sterility	0.2 µm Sterile Filtered
Packaging	HDPE Bottle
Cell Assay	BioTek Cytation 5 Imaging Plate Reader
Used as part of a Dye-Loading Buffer solution in a Cell Assay	The assay returns a passing result for the selected Dye Indicator used in the corresponding bio-validation method