

## Brilliant Sodium Assay Kit

Catalog # 9000

A total assay solution for multi-well plate based, high-throughput measurements of changes in intracellular Na<sup>+</sup> mediated through a wide variety of plasma membrane and intracellular sodium channels, sodium transporters, and non-selective cation channels.

Available in these sizes:	
9000-2	pIONeer Kit - 2 Plates
9000-10	Flex Kit - 10 Plates
9000-100	Express Kit - 100 Plates



Product Details	
Intended Use	High-throughput, no wash sodium (Na <sup>+</sup> ) assay with adherent or non-adherent cells
Compatible Instruments	Fluorescence Microscope, Plate Reader, Flow Cytometer, HTS Instrument
Product Size	Enough reagents for 2, 10, or 100 x 96-well or 384-well plate assays
Storage Temperature	Ships Ambient; Multiple Storage Conditions: -20°C and 4°C
Packaging	Components placed in bags marked with storage temperature, all packed in a mylar pouch

Kit Contents	<i>pIONeer</i> Kit Cat# 9000-2		<i>Flex</i> Kit Cat# 9000-10		<i>Express</i> Kit Cat# 9000-100	
	Size	Qty	Size	Qty	Size	Qty
ING-2 AM	50 µg Vial	2	50 µg Vial	10	500 µg Vial	10
DMSO	Not Included in pIONeer		225 µL Vial	1	2.25 mL Bottle	1
50X DySolv	800 µL Vial	1	4 mL Bottle	1	40 mL Bottle	1
10X Brilliant Assay Buffer	4 mL Bottle	1	20 mL Bottle	1	200 mL Bottle	1
50X TRS	800 µL Vial	1	4 mL Bottle	1	40 mL Bottle	1
50X Probenecid	800 µL Vial	1	4 mL Bottle	1	40 mL Bottle	1

All supplied components are tested to meet individual product specifications, including validation in a cell assay similar to that described in the No-Wash Method of the Kit Protocol, available on the product webpage. Please review the Kit Protocol for additional information and suitability of the kit with your application. Dye indicators and buffer reagents are all available for individual sale on our website. Express kits are customizable and can be tailored to your specific requirements. Contact sales for more information.