

DySolv

Catalog # 7501A

| Category | Specification |
|---|--|
| Product Details | |
| Intended Use | A solubilization aid for water insoluble dyes, used as part of a dye-loading solution |
| Product Size | 20 mL (Contact Sales for custom and bulk sized orders) |
| Product Concentration | 50X |
| Storage Temperature | 4°C |
| Product Form | |
| Appearance | Clear, Colorless solution |
| Sterility | 0.2 µm Sterile Filtered |
| Packaging | HDPE Bottle |
| Cell Assay | |
| Used as part of a Dye-Loading Buffer solution in a Cell Assay | BioTek Cytation 5 Imaging Plate Reader The assay returns a passing result for the selected Dye Indicator used in the corresponding bio-validation method |