


Fluo-4 K⁺ Salt



Membrane Impermeable form of Fluo-4, one of the most popular Green Fluorescent Calcium Indicators. Can be used in lipid membrane-free systems, liposomes, or introduced into cells via electroporation, microinjection, or other methods

Available in these sizes:

1042C	500 µg x 1
1042F	50 µg x 10



Description	
CAS number	N/A
PubChem CID	163342046
Molecular Weight	927 g/mol
Solubility	Water
K _d	355 nM
Storage	-20°C
Appearance	Light red to orange film or pellet
Container	1.5 mL polypropylene cryogenic vial

Absorbance Spectrum	Agilent Cary 60 UV-VIS Spectrophotometer
Solvent	High Calcium (Ca ²⁺) Buffer Solution
Absorbance max	494 ± 3 nm

Fluorescence Spectrum	Horiba Jobin Yvon FluoroMax Plus Spectrofluorometer
Solvent	High Calcium (Ca ²⁺) Buffer Solution
Excitation max	494 ± 3 nm
Emission max	516 ± 3 nm

LCMS	Agilent 1220 Infinity II
Column	Phenomenex 00d-4251-E0 Luna C ₁₈ , 4.6 x 100 mm, 100Å, 3 µm
UV-Vis Diode Array Detector	254 nm
Purity	≥ 90%
Single Quad MS Detector	ESI positive
Molecular Ion*	737.2 ± 0.5 m/z [M+H] ⁺ 368.6 ± 0.5 m/z [M+2H] ²⁺ *Fluo-4 K ⁺ Salt is protonated to the acid form in the mobile phase and detected in ESI Positive Mode

¹ H NMR Spectrum	Bruker Avance 400
Solvent	D ₂ O
Spectra	Product conforms to structure

In Vitro Assay	BioTek Cytation 5 Cell Imaging Multi-mode Reader
F _{Stimulus} /F _{Control} in relevant biological buffer	To Pass Test