


ING-2 TMA⁺ Salt



Description	
CAS number	N/A
PubChem CID	N/A
Molecular Weight	1087 g/mol
Solubility	Water
K _d	20 mM

LCMS	Agilent 1220 Infinity II
Column	Agilent Infinity Lab Poroshell 120 ECC18, 30 x 50 mm, 2.7 μm C ₁₈
UV-Vis Diode Array Detector	254 nm
Purity	≥ 90%
Single Quad MS Detector	ESI positive
Molecular Ion	867.27 ± 0.5 amu (MH ⁺) 889.25 ± 0.5 amu (MNa ⁺) 434.14 ± 0.5 amu (MH ₂ ²⁺)

Absorbance Spectrum	Agilent Cary 60 UV-VIS Spectrophotometer
Solvent	140 mM NaCl in MOPS
Absorbance max	517 ± 3 nm

Fluorescence Spectrum	Horiba Jobin Yvon FluoroMax 4 Spectrofluorometer
Solvent	140 mM NaCl in MOPS
Excitation max	517 nm ± 3 nm
Emission max	540 nm ± 3 nm

¹ H NMR Spectrum	Bruker Avance 400
Solvent	Deuterated water
Only relevant product peaks – with appropriate chemical shifts and peak integrations – and solvent peaks present.	